



Pride.Precision.Performance.

Improve your Sustainability Credentials AND Eliminate Lead and Lead Salts Contamination Risks by choosing Mitt-Mat®!



Mitt-Mat® outperforms conventional lead hands as the silicone acts as a supporting mechanism during flexion thus reducing metal fatigue

Mitt-Mat® offers improved value for money over conventional lead hands due to its improved durability, lower reprocessing costs and fewer replacement costs

Mitt-Mat® resistance to creasing, even in a 'doubled-over' position ensures the encapsulated lead core is not easily over stressed unlike conventional lead hands

Download the Clinical Study of
Lead Contamination on Surgical Gloves

boltons.co.uk





Pride.Precision.Performance.

- Eliminates the patient safety risks associated with lead or lead salt transfer to patients directly or via surgical gloves by completely encapsulating the lead core
- Improved flexibility means returning Mitt-Mat® to the 'flat' position for re-use is much easier
- Resistant to all commonly used cleaning chemicals
- Latex free



DESCRIPTION	CODE
MITT-MAT® ADULT	11-7070-00

What the surgeons say:

“Certainly the best of the ones that I have seen so far”

“Of all the products tested so far it has been the best”

Mitt-Mat® is available exclusively from Bolton Surgical. Please contact us or visit our website for further information.

boltons.co.uk

Bolton Surgical Limited
Churchill House, 16 Churchill Way
Chapelton, Sheffield S35 2PY
T: +44 (0)114 240 4400
E: sales@boltons.co.uk



MD 500657